

Semi Automated Dual Basket Tablet Disintegration Testing - Type PTZ AUTO 2



The PTZ AUTO Disintegration Test Instruments are available as single, double, triple and four station system which are equipped in full compliance with the USP/EP pharmacopoeia requirements. A large LCD screen and an alpha-numerical keyboard is used to edit product and test information, such as the disintegration time for each station, bath and media temperature. The total testing time for each station is displayed at the end of a run. Operator intervention to stop a run is also displayed. Print the individual disintegration time of each sample and each station when you did visually detect the disintegration of the tested sample by simply pushing the relevant character key (1-6). All this information can be printed as a test report when a PCL printer is attached to the instruments parallel printer port.

The built-in heating system uses a silent circulation pump and a heater to pre-heat the outer water bath over the range of 25-45° C; the actual bath temperature is always displayed and measured by an external manual temperature probe. The overheating protection secures automatic shut down of the pump and heater in the event that no water flow is detected. A built-in safety system including an overheating detection sensor, a thermo switch and a thermo fuse closes down the heating when temperature exceeds 60° C or max. 95°C at the heater.

For easy cleaning and maintenance, the Plexiglas water bath can be easily removed. Fast Tube Couplings connect the tubing to the heating system. The digital controlled stroke mechanism is automatically adjusted to provide the USP / EP specified 30 strokes a minute to each station. At the end of a test run, each basket is completely removed from the test media. This feature is particularly useful if sustained released products are tested, which may

PHARMA TEST AG
Siemensstrasse 5
D-63512 Hainburg (GER)

+49 6182 9532-600
+49 6182 9532-650
email@pharma-test.de
www.pharma-test.com



Semi Automated Dual Basket Tablet Disintegration Testing - Type PTZ AUTO 2

require a medium change during the run.

Upgrading

The PTZ AUTO 2 Disintegration Test Apparatus is able to be upgraded at any time to a fully automated disintegration time detection system, giving automatically independent disintegration data for each test sample using the PT-MKT auto disintegration baskets. At the end of a test run, a test report can then be printed showing the minimum and maximum disintegration time, statistics as well as bath and media temperature, stroke frequency, and product data for each sample tube and each testing station.

Featuring

- Alpha-numerical keyboard
- 8x8 cm large LC Display
- OQ/PQ Monitoring
- Test and Sample Information Filing up to 16 Methods
- Print out of Test Report at the end of a test run
- ABS - Auto-Basket-Stop of Disintegration Times, press relevant key only
- TSC - Thermostat Safety Control (includes thermo switch, thermo fuse, pressure switch, digital safety sensor)
- RS-232 and parallel Printer interface
- Built-in Basket Lift for each Testing Station
- Easy upgrading to fully automated PTZ AUTO 2EZ

Technical Data

Display:	LCD Display 8 x 8 cm
Keyboard:	alpha-numerical keys and function keys
OQ information:	programmable time sequence information to perform OQ test
Security:	password protected instrument calibration and validation procedures
Testing time:	programmable testing time from 1 second up to 59h, 59m, 59s
Test Stations:	PTZ AUTO 2 = 2 baskets, individually controllable
Thermostat:	water bath temperature programmable from 30.0 - 45.0° C, overheating protection
Accuracy:	± 0.3° C
Stroke Frequency:	30 per minute
Accuracy:	± 1 stroke
Stroke Height:	55 mm
Accuracy:	± 1 mm
Interface:	1 RS-232 serial interface for external control of the instrument and data transfer 1 parallel printer port to connect a PCL LaserJet Printer
Auto Basket Lift:	equipped with auto electric basket lift in order to raise the sample tubes out of the media at the end of a test
Upgrading:	equipped to be upgraded to a fully automated disintegration time detection system (PTZ AUTO EZ)

PHARMA TEST AG
Siemensstrasse 5
D-63512 Hainburg (GER)



+49 6182 9532-600



+49 6182 9532-650



email@pharma-test.de



www.pharma-test.com



