





UAM Motor Drive, UGD Gear and Polishing Drum PD

Polishing Drum for Sugar Coated Tablets, PD

The felt-lined polishing drum PD is designed for polishing sugar coated tablets. The covers, which lock quickly into places with three snap closures, as well as the bottom of the drums are lined with felt. All felt parts are easily exchangeable.

The cover is fitted with acrylic sheet window, which permits the observation of the polishing process. The polishing drum is suitable for screwing on to the UAM motor drive through the Universal gearing (UGD).

Technical Data:

Gross Capacity	: 2 lit
Net Capacity	: 1 lit.
Dimension	: 240 x 240 x 240 mm
Net Weight	: 6.6 Kg.
Gross Weight	: 15 Kg.
Packing	: 380 x 380 x 400 mm

We reserve the right to make technical changes without any prior notice

40-01115_PD_1.09E

Page 1

Email:email@pharmag.de Internet:www.pharmag.de



Polishing Drum PD



Universal Gear UGD

This gear is used to allow the adjustment of operating angles from 0 to about 90° and reduce the drive speed in a ratio of 1:10. Several attachments, such as Pelletisers, Cube Mixers or Coating Pans are attached via the UGD gear.

The UGD accessory has been so designed that the angle of operation can be accurately set against a scale in the instrument housing. This enables the same operational angle to be selected every time the accessory is removed and re-installed (e.g., for cleaning). Unlike some similar devices, our gear drive has been designed so that the operator can always view the on-going process as all accessories are mounted to the front of the gear drive. This allows easy viewing of the process at any time.

We reserve the right to make technical changes without any prior notice

40-01115_PD_1.09E

Email:email@pharmag.de Internet:www.pharmag.de