



Sinkers jap. PHC # 302-9000



st.st. Sinker Basket for capsules
8 mesh (approx.)
L=25x12mm # 302-9066-6
L26x9.7mm # 302-9050-6



st.st. Sinker, 4 windings, approx.
L26x9.7 mm # 302-9050-6



Sinker PTFE coated 4 windings
L26x9,7 mm # 302-9060-6
L33.5x10.4mm # 302-9062-6



Sinker PTFE coated 6 windings
L22.9x9,4 mm # 302-9061-6



st.st. Sinker 6 windings, approx.
L22.9x9.4mm # 302-9051-6
L25.4x14.6mm # 302-9053-6



stainl. steel Wire Sinker
302-9065-6



TC transdermal Cylinder incl.
Adapter (USP/EP App.6)
316-0230



40 mesh Gold pl. Basket
311-2800
10 mesh st.st. # 311-2804
40 mesh st.st. # 311-2801-6
100 mesh st.st # 311-2803-6

PHARMA TEST GmbH
Siemensstrasse 5
D-63512 Hainburg (GER)

☎ +49 6182 9532-600
☎ +49 6182 9532-650
✉ email@pharma-test.de
🌐 www.pharma-test.com



Dissolution Accessories



st.st. Stirring Shaft # 316-0116-6
st.st. Paddle Adapter # 316-0210-6



st.st. Stirring Shaft # 316-0116-6
Ctd. Paddle Adapter # 316-0250-6



st.st Basket Adapter complete
316-0220-6



Paddle over Disk - 42 mm - 68
mm - 90 mm OD (USP/EP App.5)
302-7790 (-7791 or -7792)



Ointment and Suppository
rotary Cell # 303-1015



EP Extraction Cell, perspex
303-1005



mini vessel assembly, 250 ml
303-5100-6



PTSW 0 Suppository Dialysis
Cell # 39-00502



Suppository Dissolution Basket
311-3000-6

PHARMA TEST GmbH
Siemensstrasse 5
D-63512 Hainburg (GER)



+49 6182 9532-600
+49 6182 9532-650
email@pharma-test.de
www.pharma-test.com



Dissolution Accessories



FELOPINE Baseket
303-3035



USP amber coloured vessel incl code
303-5425 (1 litre)
303-5222 (2 litre)
302-5622 (4 litre)



USP glass vessel incl. code
303-5225 (1 litre)
303-5221 (2 litre)
302-5621 (4 litre)



digital Tachometer, Handy
30-31006



digital Thermometer
30-31007



SWT 2 shaft wobble +
centricity tester # 30-31210



Height and Depth Gauge
302-9001 (paddle)
302-9002 (basket)



Depth Ball 25 mm # 315-0200



5 or 10 micron PP Sinter Filter
31-63211-50 (10 µm)
31-63511-50 (5 µm)

PHARMA TEST GmbH
Siemensstrasse 5
D-63512 Hainburg (GER)



+49 6182 9532-600
+49 6182 9532-650
email@pharma-test.de
www.pharma-test.com



Dissolution Accessories



DO² Meter # 38-31500




Manual sampling probe # 303-0100



Vibration Meter - VM
32-30511



DL1 Data-Logger for DT70 +
DT8 + PTB301/302
30-30600



Tablet Drop Magazine incl. tube
cleaning device for IDS, DTS800,
PTWS310/610/11210 # 316-6350



Spectrophotometer Calibration
Reference Cuvettes
350-2110



Flange Kit for PTFE Tubing
30-31050

We reserve the right to make technical changes without any prior notice

PHARMA TEST GmbH
Siemensstrasse 5
D-63512 Hainburg (GER)



+49 6182 9532-600
+49 6182 9532-650
email@pharma-test.de
www.pharma-test.com

